


PhD Summer School

“Biological Nitrification Inhibition: Integrating Microbial Functions, Plant Traits, and Technological Innovations for Sustainable Nitrogen Cycling”

University of Thessaly
12-16 May 2025

Organized by:



ACTIONr: Research Action Network for Reducing Reactive Nitrogen Losses from Agricultural Ecosystems

Practical Training

“Fast-track, high-throughput screening for BNIs in ammonia-oxidizing bacteria”



UNIVERSITY OF
THESSALY



This project has received funding from the European Union's Horizon 2021-2027 research and innovation programme under grant agreement No 101079299.

Contents

Introduction3

Overview of the fast-track pipeline3

Stage 2: Plate setup for sample application4

Stage 3: Plate setup for colorimetric nitrite measurements5

Data Processing.....7

Introduction

This practical session provides an introductory overview of a fast-track assay system developed to screen pure compounds and plant-derived root exudates for their biological nitrification inhibition (BNI) activity against ammonia-oxidizing microorganisms (AOM)¹. The assay includes multiple soil-relevant AOM strains; in this session, we demonstrate the use of *Nitrosospira multiformis* (ATCC25196) as a model ammonia-oxidizing bacterium (AOB).

Overview of the fast-track pipeline

Once *N. multiformis* cultures reach the late logarithmic growth phase, the fast-track assay is conducted to estimate the EC₅₀ values (Effective Concentration 50) of the test samples. The pipeline comprises three main stages (see Figure 1):

- 1) Cell harvesting and inoculation
- 2) Cell plating and sample application
- 3) Colorimetric nitrite quantification




Figure 1. Overview of the three stages of the fast-track assay.

¹ Kanellopoulos et al. 2024; BioRxiv. <https://doi.org/10.1101/2024.12.03.626636>).

In this practical training session, you will focus on stages 2 and 3 of the fast-track protocol assay. Specifically, you will: (a) prepare the “fast-track” plate, and (b) set up the plate for colorimetric nitrite quantification.

Stage 2: Plate setup for sample application

Required materials:

- One 96-well microtiter plate
- A pre-prepared 10% (v/v) concentrated culture of *N. multiformis* in Skinner-Walker (SW) medium²
- **ddH₂O control:** 1% (v/v) initial stock solution (diluted in SW medium)
- **DMSO control:** 1% (v/v) initial stock solution (diluted in SW)
- **Methyl 3-(4-Hydroxyphenyl) propionate (MHPP)** working solutions: 1% v/v DMSO initial stock solution with final concentrations in SW medium of: 0.55, 1.1, 2.75 and 5.5 mM
- **Hydrophobic root exudate working solutions:** 1% DMSO initial stock solution with final concentrations in SW medium of 0.22, 0.55, 1.1 and 1.65 µg µL⁻¹

Procedure:

1. Inoculum Dispensing

Dispense 200 µL of the 10% (v/v) *N. multiformis* concentrated culture into each well of rows B to D of the 96-well plate, as shown in Figure 2.

2. Application of samples and controls (in triplicate)

- **Negative controls:**

Add **20 µL of 1% (v/v) ddH₂O control** into wells **B1-B3** and **20 µL of 1% (v/v) DMSO control** into wells **B4-B6**.

- **Positive controls:**

Add **20 µL of MHPP working solutions** at the following final concentrations: **50, 100, 250 and 500 µM**, distributed across wells **C1-C12**.

- **Root exudate samples:**

Add **20 µL of hydrophobic root exudate working solutions** at final concentrations: **0.02, 0.05, 0.10, 0.15 µg µL⁻¹**, distributed across wells **D1-D12**.

- **Incubation**

Cover the plate with a sterile lid and incubate at **28°C for 20 hours**.

² Skinner FA, Walker N. 1961; Arch Für Mikrobiol. 38(4):339–49.

	1	2	3	4	5	6	7	8	9	10	11	12
A												
Negative controls B	ddH ₂ O	ddH ₂ O	ddH ₂ O	DMSO	DMSO	DMSO						
MHPP (positive control) C	50	50	50	100	100	100	250	250	250	500	500	500
Root exudates (samples) D	0.02	0.02	0.02	0.05	0.05	0.05	0.10	0.10	0.10	0.15	0.15	0.15
E												
F												
G												
H												

Copyright © 2009 Edita Aksamitiene

Figure 2. Plate setup for sample application.

Note:

- The final concentration of DMSO in each well will be 0.09%, a level not expected to inhibit the activity of the strain.
- Rows A and E–H are not used in this assay.

Stage 3: Plate setup for colorimetric nitrite measurements

Following the incubation period, nitrite production is quantified using the Griess reaction³, a colorimetric assay based on a two-step chemical process:

1. **Diazotizing reaction:** Nitrite reacts with **sulfanilamide** under acidic conditions to form a diazonium salt.
2. **Coupling reaction:** The diazonium salt then couples with **N-(1-naphthyl)ethylenediamine dihydrochloride (NED)** to form a **stable pink azo dye**.

The intensity of the resulting color is **directly proportional to the nitrite concentration** in the sample and is measured **spectrophotometrically at 540 nm**.

To perform this stage, you will be provided with **pre-incubated “fast-track” assay plates**, prepared identically to those assembled in Stage 2. These will be used for the colorimetric detection of nitrite accumulation, reflecting the inhibitory or stimulatory effects of the tested compounds on *N.*

³ Shinn MB. ACS Publications. American Chemical Society; 1941. <https://pubs.acs.org/doi/abs/10.1021/i560089a010>

multiformis activity.

Required materials:

- **Incubated 96-well plate** from Stage 2 (containing the test samples)
- **One fresh 96-well microtiter plate** (for colorimetric measurement)
- **Sodium nitrite (NaNO_2) standard solutions** at the following concentrations: 0, 1, 2, 5, 10, 20, 50, 100 and 200 μM
- **Skinner and Walker (SW) medium** in a **30 mL sterile polystyrene Sterilin™ bottle**
- **Diazotizing reagent** (typically contains sulfanilamide in acidic solution)
- **Coupling reagent** (typically contains N-(1-naphthyl)ethylenediamine dihydrochloride, NED)

Procedure:

1. Preparation of standard curve

Add 100 μL of NaNO_2 standard curve solutions in the range: 0, 1, 2, 5, 10, 20, 50, 100 and 200 μM into wells A1 to A9 of the 96-well plate to establish the nitrite standard curve.

2. Sample dilution

Dispense 90 μL of SW medium into the appropriate wells of rows B to D (see Figure 3). Then, transfer 10 μL of each sample from the pre-incubated assay plate (from stage 2) into the corresponding wells of the new plate (rows B-D). This results in a 1:10 dilution of each sample to ensure compatibility with the standard curve range.

3. Color development via Griess reaction

Add **20 μL of diazotizing reagent** to each well, followed by **20 μL of coupling reagent**. A magenta color will develop immediately in wells containing nitrite, with **color intensity directly proportional to the nitrite concentration**.

4. Incubation

Incubate the plate for **10 minutes at room temperature** to allow full color development.

5. Nitrite measurement

Measure the **absorbance at 540 nm** using a microplate reader.

	1	2	3	4	5	6	7	8	9	10	11	12
NaNO ₂ standard solutions	A	0	1	2	5	10	20	50	100	200		
Negative controls	B	ddH ₂ O	ddH ₂ O	ddH ₂ O	DMSO	DMSO	DMSO					
MHPP (positive control)	C	50	50	50	100	100	100	250	250	250	500	500
Root exudates (samples)	D	0.02	0.02	0.02	0.05	0.05	0.05	0.10	0.10	0.10	0.15	0.15
Dilution 1:10												
	F											
	G											
	H											

Copyright © 2009 Edita Aksamitiene

Figure 3. Plate setup for the nitrite measurements

Data Processing

1. Normalize nitrite values

For each sample and control, normalize the nitrite concentration by dividing by the average value of the negative controls:

$$\text{Normalised } NO_{2\text{control}} = \frac{NO_{2\text{control}}}{\text{Average } NO_{2\text{control}}}$$

$$\text{Normalised } NO_{2\text{sample}} = \frac{NO_{2\text{sample}}}{\text{Average } NO_{2\text{control}}}$$

2. Calculate Ammonia-Oxidation Inhibition (AOI%)

For each replicate, calculate the AOI% based on the deviation from the control, following the equation:

$$AOI \% = - (\text{Normalised } NO_{2\text{sample}} - \text{Normalised } NO_{2\text{control}})$$

3. Plot the AOI% vs. Concentration

Use the AOI% values for each tested concentration to generate inhibition curves. These plots can be used to estimate the EC₅₀ (Effective Concentration at which 50% inhibition is observed), using curve-fitting methods (e.g., non-linear regression).